

Product Datasheet

LRRC46 antibody (orb1124483)

Description	Rabbit polyclonal antibody to LRRC46
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Rabbit polyclonal LRRC46 (1) antibody was raised against a recombinated human LRRC46 protein 4-166aa (BC010412).
Preservatives	PBS, pH 7.4 with 0.02% Sodium Azide
Form/Appearance	Liquid
Concentration	100ug/100ul
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.