

Product Datasheet

ZNF10 antibody (orb1124611)

Description	Rabbit polyclonal antibody to ZNF10
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Rabbit polyclonal ZNF10 (1) antibody was raised against a recombinant human ZNF10 protein 20-242aa (BC024182)..
Preservatives	PBS, pH 7.4 with 0.02% Sodium Azide
Form/Appearance	Liquid
Concentration	100ug/100ul
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.