

Product Datasheet

c-Myc antibody (orb1125074)

Description	Mouse monoclonal antibody to c-Myc
Species/Host	Mouse
Reactivity	Other
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Synthetic peptide corresponding to fragment (aa/408-420) c-Myc
Preservatives	PBS, pH 7.4 with 0.02% sodium azide.
Form/Appearance	Liquid
Concentration	100ug/100ul
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Application notes	"WB (1:1000) analysis of c-Myc tagged rice protein A in leaf and panicle separated by 10% SDS-PAGE with Anti-c-Myc (AbM59008-2B2-PU)."
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	"WB (1:1000) analysis of c-Myc tagged rice protein A in leaf and panicle separated by 10% SDS-PAGE with Anti-c-Myc
Expiration Date	12 months from date of receipt.