

Product Datasheet

MKI67 antibody (orb1125122)

Description	Mouse monoclonal antibody to MKI67
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Recombinant protein corresponding to full length MKI67.
Preservatives	PBS, pH 7.4 with 0.02% sodium azide.
Form/Appearance	Liquid
Concentration	100ug/100ul
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Application notes	IHC (1:5000, 20X) analysis of MKI67 expression in breast cancer tissue with Anti-MKI67 (AbM51101-2-PU).
Isotype	IgG2a
Clonality	Monoclonal
Dilution Range	IHC (1:5000, 20X) analysis of MKI67 expression in breast cancer tissue with Anti-MKI67
Expiration Date	12 months from date of receipt.