

Product Datasheet

HSPA8 antibody (orb1125317)

Description	Mouse monoclonal antibody to HSPA8
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant protein corresponding to full length HSPA8.
Preservatives	citrate-tris-HCl buffer pH 7.0 with 0.02% Proclin 303.
Form/Appearance	Liquid
Concentration	100ug/100ul
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Application notes	None
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	None
Expiration Date	12 months from date of receipt.