

## Product Datasheet

### Uncoupling Protein 1 (UCP1) antibody [Out of stock] (orb1125416)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | Rabbit polyclonal antibody to Uncoupling Protein 1 (UCP1)  |
| <b>Species/Host</b>        | Rabbit   |
| <b>Reactivity</b>          | Human, Mouse, Other, Rat   |
| <b>Conjugation</b>         | Unconjugated   |
| <b>Tested Applications</b> | ELISA, IHC-P, WB   |
| <b>Immunogen</b>           | 12-aa peptide from Mouse/ Rat UCP1 - KLH   |
| <b>Form/Appearance</b>     | Liquid   |
| <b>Storage</b>             | This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt. |
| <b>Note</b>                | For research use only  |
| <b>Clonality</b>           | Polyclonal   |
| <b>Purity</b>              | Affinity purified  |
| <b>Dilution Range</b>      | Western Blot 1 -10µg/ml. ELISA 1:10,000 to 1: 100,000. IHC at 2-20µg/ml in formaldehyde fixed tissue   |
| <b>Expiration Date</b>     | 12 months from date of receipt.  |