

Product Datasheet

Poly-histidine Tag (6-His) antibody [Out of stock] (orb1125442)

Description	Rabbit polyclonal antibody to Poly-histidine Tag (6-His)
Species/Host	Rabbit
Reactivity	Other
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	6-His (HHHHHH) conjugated to KLH
Preservatives	PBS, pH 7.4 with 0.02% sodium azide.
Form/Appearance	Liquid
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purified
Dilution Range	ELISA: Coating 1:100-1:500, Primary 1:1,000-1:30,000. ImmunoCytochemistry: 1:100-1:400. Western Blot 1:1,000 - 1:30,000
Expiration Date	12 months from date of receipt.