

Product Datasheet

Green Fluorescent Protein, (GFP) antibody [Out of stock] (orb1125519)

Description

Goat polyclonal antibody to Green

Species/Host

Goat

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IHC-P, WB

Applications
Immunogen

GST-GFP fusion protein corresponding to the full length amino acid sequence (246aa) derived from the jellyfish Aequorea victoria.

Preservatives

Liquid sterile filtered; 0.02M Potassium Phosphate, 0.15M Sodium Chloride, pH7.2 and 0.01% Sodium Azide.

Form/Appearance

Liquid

Storage

This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

Note

For research use only

Clonality

Polyclonal

Purity

Affinity purified

Dilution Range

Suitable for ELISA (1:4,000-1:16,000), immunoblotting (Western 1:400-1:2,000), and immunohistochemistry (1:200-1:1,000).

Expiration Date

12 months from date of receipt.