

Product Datasheet

**Estrogen Receptor beta antibody [Out of stock]
(orb1125525)**

Description

Rabbit polyclonal antibody to Estrogen

Species/Host

Rabbit

Reactivity

Human, Other

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications
Immunogen

Synthetic peptide corresponding to the n-terminal amino acid residues 55-70 of rat and Mouse ERB conjugates to KLH. (A55 EPQKSPWCEARSLEH70)

Preservatives

50 µg of affinity Purified IgG in 50 µl PBS with 1 mg/ml rAlbumin and 0.05% Sodium Azide

Form/Appearance

Liquid

Storage

This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

Note

For research use only

Clonality

Polyclonal

Purity

Affinity purified

Dilution Range

Suggested starting dilutions :
Immunohistochemistry 1-2 µg/ml Gel
Supershift: 1-5 µl Western Blot 1-2 µg/ml (band approx 55kDa)

Expiration Date

12 months from date of receipt.