

Product Datasheet

Green Fluorescent Protein, (GFP) antibody [Out of stock] (orb1126232)

Description

Rabbit polyclonal antibody to Green

Species/Host	Rabbit
Reactivity	Other
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	GST-GFP fusion protein corresponding to the full length amino acid sequence (246aa) derived from jellyfish Aequorea victoria.
Preservatives	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride , pH 7.2 with 0.01% Gentamicin Sulfate.
Form/Appearance	Liquid
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purified
Dilution Range	Recommended Dilutions 1:10 000-1:50 000 (ELISA), 1:1000-1:5000 (Western Blot), 1:500-1:3000 (Immunohistochemistry)
Expiration Date	12 months from date of receipt.