

Product Datasheet

Helicobacter pylori antibody [Out of stock] (orb1126270)

Description	Rabbit polyclonal antibody to Helicobacter pylori
Species/Host	Rabbit
Reactivity	Other
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Helicobacter pylori
Preservatives	PBS buffer with 1% rAlbumin and Thimerosal
Form/Appearance	Liquid
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purified
Dilution Range	Recommended Dilution 1:1000 (EIA), 1:200 (Immunoblot)
Expiration Date	12 months from date of receipt.