



Product Datasheet

MDR1 antibody (orb11267)


Description	MDR1 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human MDR1. Please contact us for the exact immunogen sequence. The peptide is available as orb374742.
Target	MDR1
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	141 kDa
Uniprot ID	P08183, P06795, P43245
NCBI	21, 000927, 000918, 41
Entrez	5243
Dilution Range	IF/ICC: 1:1000, IHC-P:1:1000 (based on 0.5 mg/ml)
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...



Immunohistochemical analysis of formalin...



IHC-P of human lung cancer tissue(p-Glyc...