

Product Datasheet

Fibroblast Activation Protein alpha (FAP) antibody [Out of stock] (orb1127121)

Description	Rabbit polyclonal antibody to Fibroblast Activation Protein alpha (FAP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Synthetic peptide - KLH conjugated
Form/Appearance	Liquid
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purified
Expiration Date	12 months from date of receipt.