

Product Datasheet

Prostatic Acid Phosphatase (PAP) antibody [Out of stock] (orb1127275)

Description

Rabbit polyclonal antibody to

Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	RIA
Immunogen	Affinity Purified PAP
Preservatives	Neat serum in 0.1% sodium azide.
Form/Appearance	Liquid
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Isotype	Neat Serum
Clonality	Polyclonal
Purity	Affinity purified
Dilution Range	We recommend that you optimize titer using the reagents in your unique assay system. 70,000 tt/ml in RIA DAB. Lot to lot consistency of a polyclonal pool is about 25%
Expiration Date	12 months from date of receipt.