

Product Datasheet

Phosphorylated Tau (PHF-Tau) antibody [Out of stock] (orb1127877)

Description

Mouse monoclonal antibody to

Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	human PHF-Tau
Preservatives	no carrier proteins or preservatives added
Form/Appearance	Liquid
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Affinity purified
Dilution Range	Suggested use: Western Blot @ 20-60µg/ml will detect 50ngSDS-denatured and bmetcaptoethanol reduced PHF-Tau with chemiluminescence detection. ELISA @ 5-10µg/ml as capture reagent. Histochemistry @ 5-10µg/ml in formalin fixed paraffin embedded brain tissue. This product must be titrated in end-users assay. Endotoxin <0.0025 EU/µg
Expiration Date	12 months from date of receipt.