

Product Datasheet

Synaptopodin (Podocyte Protein) antibody [Out of stock] (orb1127881)

Description	Mouse monoclonal antibody to Synaptopodin (Podocyte Protein)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, WB
Immunogen	Isolated rat Kidney glomeruli
Preservatives	PBS, pH 7.4 with 0.02% sodium azide.
Form/Appearance	Liquid
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Affinity purified
Dilution Range	Immunofluorescence and immunoblotting (approx 1:10-1:50)
Expiration Date	12 months from date of receipt.