

## Product Datasheet

### Polyhistidine Tag (HIS6-Tag) antibody [Out of stock] (orb1127927)

|                            |                                                                                                                                                                                                                                                                                |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody to Polyhistidine Tag (HIS6-Tag)                                                                                                                                                                                                                      |
| <b>Species/Host</b>        | Mouse                                                                                                                                                                                                                                                                          |
| <b>Reactivity</b>          | Other                                                                                                                                                                                                                                                                          |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                                                                                                                                                   |
| <b>Tested Applications</b> | ELISA, IP, WB                                                                                                                                                                                                                                                                  |
| <b>Preservatives</b>       | PBS, pH 7.4 with 0.02% sodium azide.                                                                                                                                                                                                                                           |
| <b>Form/Appearance</b>     | Liquid                                                                                                                                                                                                                                                                         |
| <b>Storage</b>             | This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt. |
| <b>Note</b>                | For research use only                                                                                                                                                                                                                                                          |
| <b>Clonality</b>           | Monoclonal                                                                                                                                                                                                                                                                     |
| <b>Purity</b>              | Affinity purified                                                                                                                                                                                                                                                              |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                                                                                                                                                |