

Product Datasheet

NA⁺/K⁺ ATPase alpha-1 antibody [Out of stock] (orb1127954)

Description	Mouse monoclonal antibody to NA ⁺ /K ⁺ ATPase alpha-1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Immunogen	Purified NA ⁺ / K ⁺ ATPase isolated from membrane fractions of rabbit kidney outer medulla
Preservatives	0.1M Tris Glycine, pH 7.4, 0.15 M NaCl and 0.05% Sodium Azide
Form/Appearance	Liquid
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Clonality	Monoclonal
Purity	Affinity purified
Expiration Date	12 months from date of receipt.