

Product Datasheet

Green Fluorescent Protein, (GFP) antibody [Out of stock] (orb1128085)

Description

Mouse monoclonal antibody to

Species/Host	Mouse
Reactivity	Other
Conjugation	Unconjugated
Tested Applications	Cell Culture, ELISA, IF, IP, WB
Immunogen	Full length amino acid sequence (246aa) derived from the jelly fish <i>Aequorea victoria</i>
Preservatives	Prepared from tissue culture supernatant containing 0.02M Sodium Phosphate, 0.15M Sodium Chloride; pH 7.2
Form/Appearance	Liquid
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Isotype	γ 2a
Clonality	Monoclonal
Purity	Affinity purified
Dilution Range	ICC/IF, IP, WB, ELISA
Expiration Date	12 months from date of receipt.