

Product Datasheet

**Adeno Associated Virus (AAV), Rep Protein
antibody [Out of stock] (orb1128458)**

Description

Mouse monoclonal antibody to

	Rep Protein
Species/Host	Mouse
Reactivity	Human, Other
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant AAV-2 Rep 78 protein, N-terminally truncated by 171 aa.
Preservatives	50µg (lyophilized); contains 0.09% sodium azide and 0.5% rAlbumin.
Form/Appearance	Liquid
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Affinity purified
Dilution Range	Suitable for immunoblotting. Working dilution: reconstitute with 1ml distilled water and dilute 1:50-100 for immunoblotting. Dilution buffer: PBS or Tris-buffered saline.
Expiration Date	12 months from date of receipt.