

Product Datasheet

Adeno Associated Virus (AAV) replicase antibody
[Out of stock] (orb1128464)

Description

Mouse monoclonal antibody to

replicase

Species/Host

Mouse

Reactivity

Human, Other

Conjugation

Unconjugated

Tested

IF, IP

Applications
Immunogen

Recombinant AAV-2 Rep78 protein, N-terminally truncated by 171 aa

Preservatives

Contains 0.09% sodium azide and 0.5% rAlbumin. Reconstitute in 1 ml of distilled water. For further dilutions use PBS or Tris.

Form/Appearance

Liquid

Storage

This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

Note

For research use only

Isotype

IgG1

Clonality

Monoclonal

Purity

SpASS Protein A recombinant.

Dilution Range

Suitable for immunofluorescence (1:10) and immunoprecipitation. Working Dilution: After reconstitution with 1ml distilled water, dilute 1:10 (immunofluorescence).

Expiration Date

12 months from date of receipt.