

Product Datasheet

Prostatic Acid Phosphatase (PAP) antibody [Out of stock] (orb1128494)

Description

Mouse monoclonal antibody to

Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	Cell Culture, IF, IHC-P
Immunogen	Highly pure human PAP
Preservatives	supplied in PBS buffer, pH 7.5 with 0.1% NaN ₃
Form/Appearance	Liquid
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified (mouse ascites) by Protein A affinity chromatography
Dilution Range	ICC/IF: Use at an assay dependent dilution. IHC-P: 1/200. Do not perform antigen retrieval.
Expiration Date	12 months from date of receipt.