

Product Datasheet

Helicobacter pylori antibody [Out of stock] (orb1128636)

Description	Mouse monoclonal antibody to Helicobacter pylori
Species/Host	Mouse
Reactivity	Other
Conjugation	Unconjugated
Immunogen	Purified H. pylori antigen
Preservatives	PBS buffer, pH 7.0, with 0.02% sodium azide.
Form/Appearance	Liquid
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Affinity purified
Expiration Date	12 months from date of receipt.