

Product Datasheet

Ras antibody (orb1129165)

Description	Rat monoclonal antibody to Ras
Species/Host	Rat
Reactivity	Human, Mouse, Other, Rat
Conjugation	Unconjugated
Tested Applications	IP
Immunogen	Syngeneic normal rat kidney cells (NRK) that had been infected with Harvey murine sarcoma virus.
Preservatives	PBS, pH 7.4 with 0.02% sodium azide.
Form/Appearance	Liquid
Concentration	100ug/100ul
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.