

Product Datasheet

RICK antibody (orb1129264)

Description

Rabbit polyclonal antibody to RICK

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

WB

Applications
Immunogen

Synthetic peptide derived from Human human RICK.

Preservatives

PBS, pH 7.4 with 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1ug/ul

Storage

This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

Note

For research use only

Application notes

Recommended dilutions WB: 1/2000. Detects a band of approximately 65 kDa (predicted molecular weight: 61 kDa). Not yet tested in other applications. Optimal dilutions/concentrations should be determined by the end user.

Isotype

IgG

Clonality

Polyclonal

Dilution Range

Recommended dilutions WB: 1/2000. Detects a band of approximately 65 kDa (predicted molecular weight: 61 kDa). Not yet tested in other applications. Optimal dilutions/concentrations should be determined by the end user.

Expiration Date

12 months from date of receipt.