

Product Datasheet

Sortilin (neurotensin receptor 3) antibody [Out of stock] (orb1129459)

Description

Rabbit polyclonal antibody to

Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P
Immunogen	Extracellular domain of glycosylated human sortilin produced in CHO cells was used as the immunogen.
Preservatives	PBS, pH 7.4 with 0.02% sodium azide.
Form/Appearance	Liquid
Concentration	1ug/ul
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Clonality	Polyclonal
Purity	Filter purified
Dilution Range	A dilution of 1:1000 is recommended for IHC. The optimal working dilution should be determined by the end user.
Expiration Date	12 months from date of receipt.