

Product Datasheet

capsaicin receptor (608-621) antibody [Out of stock] (orb1129497)

Description

Rabbit polyclonal antibody to

Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A synthetic peptide (PSESTSHRWGPA) of human capsaicin receptor protein (aa: 608-621) conjugated to KLH has been used as the antigen.
Preservatives	PBS, pH 7.4 with 0.02% sodium azide.
Form/Appearance	Liquid
Concentration	1ug/ul
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Clonality	Polyclonal
Purity	Filter purified
Dilution Range	Use at 1:1000 to 1:2000 dilution. This antibody has not been tested in other applications. The optimal dilution should be determined by the end user.
Expiration Date	12 months from date of receipt.