

Product Datasheet

Capsaicin Receptor (531-541) antibody [Out of stock] (orb1129506)

Description

Rabbit polyclonal antibody to

Species/Host

Rabbit

Reactivity

Mouse, Rat

Conjugation

Unconjugated

Tested

IHC-P

Applications
Immunogen

A synthetic peptide (YFSLKEYVAS) of human capsaicin receptor protein (aa: 531-541) conjugated to KLH has been used as the immunogen.

Preservatives

PBS, pH 7.4 with 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1ug/ul

Storage

This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

Note

For research use only

Clonality

Polyclonal

Purity

Filter purified

Dilution Range

Use at 1:1000 to 1:2000 dilution. This antibody has not been tested in other applications. The optimal dilution should be determined by the end user.

Expiration Date

12 months from date of receipt.