

Product Datasheet

rat capsaicin receptor (819-838) antibody [Out of stock] (orb1129510)

Description

Mouse monoclonal antibody to rat

Species/Host	Mouse
Reactivity	Mouse, Other, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	A synthetic peptide (C-GSLKPEDAEVFKDSMPGK) as a part of rat VR1 protein (aa: 819-838) conjugated to KLH has been used as the immunogen.
Preservatives	PBS, pH 7.4 with 0.02% sodium azide.
Form/Appearance	Liquid
Concentration	1ug/ul
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	Protein G Purified
Dilution Range	Recommended to be used at a concentration of 0.5-1 µg/ml. The optimal dilution should be determined by the end user.
Expiration Date	12 months from date of receipt.