

Product Datasheet

Potassium Channel KCa2.3, (C-term) (Small-conductance Ca²⁺-activated Potassium Channel 3) antibody [Out of stock] (orb1130077)

Description

Rabbit polyclonal antibody to

(Small-conductance Ca²⁺-activated P...

Species/Host

Rabbit

Reactivity

Human, Mouse, Other, Rat

Conjugation

Unconjugated

Tested

WB

Applications
Immunogen

Synthetic peptide corresponding to amino acid residues 659-674 of human KCa2.3. This sequence has 100% homology in rat, mouse, pig.

Preservatives

PBS, pH 7.4 with 0.05% sodium azide.

Form/Appearance

Liquid

Concentration

1ug/ul

Storage

This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

Note

For research use only

Clonality

Polyclonal

Purity

affinity isolated antibody

Expiration Date

12 months from date of receipt.