

Product Datasheet

Glutamate Receptor 4, Metabotropic antibody [Out of stock] (orb1130165)

Description

Rabbit polyclonal antibody to

Metabotropic

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P

Applications
Immunogen

synthetic peptide corresponding to the N-terminal extracellular domain of human metabotropic glutamate receptor 4, conjugated to KLH. The immunizing peptide has 100% homology with the rat gene and 89% homology with the mouse gene.

Preservatives

PBS, pH 7.4 with 0.05% sodium azide.

Form/Appearance

Liquid

Concentration

1ug/ul

Storage

This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

Note

For research use only

Clonality

Polyclonal

Purity

affinity isolated antibody

Expiration Date

12 months from date of receipt.