

## Product Datasheet

### phospho-GSK-3 beta (pSer9) antibody [Out of stock] (orb1130507)

<b>Description</b>	Rabbit polyclonal antibody to phospho-GSK-3 beta (pSer9)
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, WB
<b>Immunogen</b>	synthetic phosphopeptide derived from the region of GSK-3 $\beta$ that contains Ser9.
<b>Preservatives</b>	PBS, pH 7.4 with 0.05% sodium azide.
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	1 $\mu$ g/ $\mu$ l
<b>Storage</b>	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Purity</b>	affinity isolated antibody
<b>Expiration Date</b>	12 months from date of receipt.