

Product Datasheet

MEK1 (phospho-Thr292) from antibody [Out of stock] (orb1130532)

Description	Rabbit polyclonal antibody to MEK1 (phospho-Thr292) from
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	synthetic phosphopeptide derived from the region of MEK1 that contains threonine 292.
Preservatives	PBS, pH 7.4 with 0.05% sodium azide.
Form/Appearance	Liquid
Concentration	1ug/ul
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Clonality	Polyclonal
Purity	affinity isolated antibody
Expiration Date	12 months from date of receipt.