

Product Datasheet

phospho-S6-Kinase (pThr421/pSer424) antibody [Out of stock] (orb1130568)

Description	Rabbit polyclonal antibody to phospho-S6-Kinase (pThr421/pSer424)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	synthetic double phospho-Thr421/Ser424 peptide corresponding to residues around Thr421/Ser424 of human p70S6K.
Preservatives	PBS, pH 7.4 with 0.05% sodium azide.
Form/Appearance	Liquid
Concentration	1ug/ul
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Clonality	Polyclonal
Purity	affinity isolated antibody
Expiration Date	12 months from date of receipt.