

Product Datasheet

beta -Actin-Peroxidase antibody (orb1131388)

Description	Mouse monoclonal antibody to beta - Actin-Peroxidase
Species/Host	Mouse
Reactivity	Aves, Bovine, Human, Mouse, Other, Rat, Sheep
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	slightly modified β -cytoplasmic actin N-terminal peptide, Ac-Asp-Asp-Asp-Ile-Ala-Ala-Leu-Val-Ile-Asp-Asn-Gly-Ser-Gly-Lys, conjugated to KLH.
Preservatives	PBS, pH 7.4 with 0.05% sodium azide.
Form/Appearance	Liquid
Concentration	1ug/ul
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at - 20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.