

Product Datasheet

Tubulin, Acetylated antibody [Out of stock] (orb1131534)

Description	Mouse monoclonal antibody to Tubulin, Acetylated
Species/Host	Mouse
Reactivity	Bovine, Human, Mouse, Other, Rat, Virus
Conjugation	Unconjugated
Tested Applications	IHC-P, RIA, WB
Immunogen	acetylated β -tubulin from the outer arm of sea urchin sperm axonemes.
Preservatives	PBS, pH 7.4 with 0.05% sodium azide.
Form/Appearance	Liquid
Concentration	1ug/ul
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Purity	ascites fluid
Expiration Date	12 months from date of receipt.