

Product Datasheet

RanGAP1 (C-terminal) antibody [Out of stock] (orb1131741)

Description	Rabbit polyclonal antibody to RanGAP1 (C-terminal)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	synthetic peptide corresponding to amino acids 573-587 located at the C-terminus of human RanGAP1. The sequence is highly conserved (80% identity) in mouse RanGAP1.
Preservatives	PBS, pH 7.4 with 0.05% sodium azide.
Form/Appearance	Liquid
Concentration	1ug/ul
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Clonality	Polyclonal
Purity	affinity isolated antibody
Expiration Date	12 months from date of receipt.