

Product Datasheet

phospho-eIF4G (pSer1108) antibody [Out of stock] (orb1131811)

| | |
|----------------------------|--|
| Description | Rabbit polyclonal antibody to phospho-eIF4G (pSer1108) |
| Species/Host | Rabbit |
| Reactivity | Human, Other |
| Conjugation | Unconjugated |
| Tested Applications | WB |
| Immunogen | synthetic phosphopeptide derived from the region of human eIF4G that contains serine 1108. Rabbit is 92% homologous. |
| Preservatives | PBS, pH 7.4 with 0.05% sodium azide. |
| Form/Appearance | Liquid |
| Concentration | 1ug/ul |
| Storage | This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt. |
| Note | For research use only |
| Clonality | Polyclonal |
| Purity | affinity isolated antibody |
| Expiration Date | 12 months from date of receipt. |