

Product Datasheet

Glutamic Acid Decarboxylase 65 (5-22) antibody
[Out of stock] (orb1132108)

Description

Rabbit polyclonal antibody to

(5-22)

Species/Host

Rabbit

Reactivity

Human, Mouse, Other, Rat

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications
Immunogen

synthetic peptide corresponding to the N-terminal region of human GAD 65 (amino acids 5-22). The sequence is highly conserved in rat and pig GAD 65 (single amino acid substitution) and mouse and monkey GAD 65, but is not found in GAD 67.

Preservatives

PBS, pH 7.4 with 0.05% sodium azide.

Form/Appearance

Liquid

Concentration

1ug/ul

Storage

This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

Note

For research use only

Clonality

Polyclonal

Purity

IgG fraction of antiserum

Expiration Date

12 months from date of receipt.