

Product Datasheet

Glutamic Acid Decarboxylase 65/67 antibody [Out of stock] (orb1132111)

Description

Rabbit polyclonal antibody to

65/67

Species/Host

Rabbit

Reactivity

Human, Mouse, Other, Rat

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications
Immunogen

synthetic peptide corresponding to the C-terminal region of human GAD 67 (amino acids 579-594). The sequence is identical in human GAD 65 (amino acids 570-585), in rat, mouse, and pig GAD 65, in rat, pig, and cat GAD 67, and is highly conserved in GAD 67 of mouse origin (single amino acid substitution).

Preservatives

PBS, pH 7.4 with 0.05% sodium azide.

Form/Appearance

Liquid

Concentration

1ug/ul

Storage

This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

Note

For research use only

Clonality

Polyclonal

Purity

IgG fraction of antiserum

Expiration Date

12 months from date of receipt.