

Product Datasheet

Glutamine Synthetase antibody [Out of stock] (orb1132112)

Description	Rabbit polyclonal antibody to Glutamine Synthetase
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	synthetic peptide corresponding to the C-terminal region of mouse glutamine synthetase (GS) (amino acids 357-373. This sequence is identical in human, bovine, rat, hamster, and pig GS.
Preservatives	PBS, pH 7.4 with 0.05% sodium azide.
Form/Appearance	Liquid
Concentration	1ug/ul
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Clonality	Polyclonal
Purity	IgG fraction of antiserum
Expiration Date	12 months from date of receipt.