

Product Datasheet

Neurofilament-M O-glycosylated antibody [Out of stock] (orb1132149)

Description

Mouse monoclonal antibody to

Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	Cell Culture, IF, IHC-P, WB
Immunogen	protein fraction enriched for components of the human neuronal cytoskeleton from cell lines NT2-N and SK-N-BE(2). The epitope is localized in the carboxy-terminal tail domain of NF-M.
Preservatives	PBS, pH 7.4 with 0.05% sodium azide.
Form/Appearance	Liquid
Concentration	1ug/ul
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Clonality	Monoclonal
Purity	purified immunoglobulin
Expiration Date	12 months from date of receipt.