

Product Datasheet

Synaptopodin (RS-18) antibody [Out of stock] (orb1132175)

Description	Rabbit polyclonal antibody to Synaptopodin (RS-18)
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	synthetic peptide corresponding to amino acids 431-448 of rat synaptopodin conjugated to KLH. This rat synaptopodin sequence has ~80% identity with human synaptopodin.
Preservatives	PBS, pH 7.4 with 0.05% sodium azide.
Form/Appearance	Liquid
Concentration	1ug/ul
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Clonality	Polyclonal
Purity	IgG fraction of antiserum
Expiration Date	12 months from date of receipt.