

Product Datasheet

Vesicular Acetylcholine Transporter (VACHT) antibody (orb1132189)

Description	Rabbit polyclonal antibody to Vesicular Acetylcholine Transporter (VACHT)
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	synthetic peptide corresponding to the C-terminal of rat vesicular acetylcholine transporter (VACHT) (amino acids 512-530).
Preservatives	PBS, pH 7.4 with 0.05% sodium azide.
Form/Appearance	Liquid
Concentration	1ug/ul
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.