

Product Datasheet

Oligodendrocyte Marker O4 antibody [Out of stock] (orb1132378)

Description	Mouse monoclonal antibody to Oligodendrocyte Marker O4
Species/Host	Mouse
Reactivity	Human, Mouse, Rat, Virus
Conjugation	Unconjugated
Tested Applications	FC, IHC-P
Immunogen	white matter of corpus callosum from bovine brain.
Preservatives	PBS, pH 7.4 with 0.05% sodium azide.
Form/Appearance	Liquid
Concentration	1ug/ul
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Isotype	IgM
Clonality	Monoclonal
Purity	purified immunoglobulin
Expiration Date	12 months from date of receipt.