

Product Datasheet

**KCa2.3 (C-term) (SK3) antibody [Out of stock]
(orb1132413)**

Description

Rabbit polyclonal antibody to KCa2.3

Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Peptide (C)KQIGSLESKLEHTAS, corresponding to amino acid residues 659-674 of human KCa2.3 (Accession Q9UGI6).
Preservatives	PBS, pH 7.4 with 0.05% sodium azide.
Form/Appearance	Liquid
Concentration	1ug/ul
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purified on immobilized antigen
Dilution Range	Western Blotting: Rat brain membranes (1:200)
Expiration Date	12 months from date of receipt.