



## Product Datasheet

### c-Raf (phospho-Ser338/Tyr340) antibody (orb11325)


|                            |                                                                                                                                      |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | c-Raf (phospho-Ser338/Tyr340) antibody                                                                                               |
| <b>Species/Host</b>        | Rabbit                                                                                                                               |
| <b>Reactivity</b>          | Human, Mouse                                                                                                                         |
| <b>Conjugation</b>         | Unconjugated                                                                                                                         |
| <b>Tested Applications</b> | FC, WB                                                                                                                               |
| <b>Immunogen</b>           | KLH conjugated Synthesised phosphopeptide derived from human RAF1 around the phosphorylation site of Ser338/Tyr340 (RD(p-S)S(p-Y)YW) |
| <b>Target</b>              | phospho-RAF1 (Ser338/Tyr340)                                                                                                         |
| <b>Preservatives</b>       | 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.                                                                |
| <b>Form/Appearance</b>     | Liquid                                                                                                                               |
| <b>Concentration</b>       | 1mg/ml                                                                                                                               |
| <b>Storage</b>             | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.                                                     |
| <b>Note</b>                | For research use only                                                                                                                |
| <b>Isotype</b>             | IgG                                                                                                                                  |
| <b>Clonality</b>           | Polyclonal                                                                                                                           |
| <b>MW</b>                  | 73 kDa                                                                                                                               |
| <b>Uniprot ID</b>          | <b>P04049</b>                                                                                                                        |
| <b>Dilution Range</b>      | WB=1:500-2000, Flow-Cyt=2ug/Test                                                                                                     |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                      |



Immunohistochemical staining of human glomus tissue.



Western blot analysis of mouse heart lysate.



rat brain tissue (Lane 2) using c-Raf (p-Ser338/Tyr340) antibody.