

Product Datasheet

**Cav2.3 (alpha 1E, R-type) antibody [Out of stock]
(orb1132507)**

Description

Rabbit polyclonal antibody to Cav2.3

Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Peptide (C)SASQ ERSLD EGVSI DG (CNE2), corresponding to residues 892-907 of rat CaV2.3 (Accession Q07652).
Preservatives	PBS, pH 7.4 with 0.05% sodium azide.
Form/Appearance	Liquid
Concentration	1ug/ul
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Clonality	Polyclonal
Purity	affinity purified on immobilized antigen.
Dilution Range	Western Blotting: Rat brain membranes (1:200). Immunohistochemistry: Rat brain sections.
Expiration Date	12 months from date of receipt.