

Product Datasheet

glycated hemoglobin HbA1c antibody [Out of stock] (orb1132543)

Description	Mouse monoclonal antibody to glycated hemoglobin HbA1c
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Immunogen	Glycated (in vitro) synthetic peptide corresponding to the N-terminal end of beta-globin.
Preservatives	PBS, pH 7.4 with 0.05% sodium azide.
Form/Appearance	Liquid
Concentration	0.1 mg/ml
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Isotype	IgG for MAb
Clonality	Monoclonal
Purity	Size-exclusion chromatography
Expiration Date	12 months from date of receipt.