

Product Datasheet

Helicobacter pylori CagA-protein antibody [Out of stock] (orb1132578)

Description	Mouse monoclonal antibody to Helicobacter pylori CagA-protein
Species/Host	Mouse
Reactivity	Other
Conjugation	Unconjugated
Immunogen	Recombinant fragment of CA-gene of H. pylori consisting of 580 C-terminal amino acid residues (calculated MW ~ 70 kDa)
Preservatives	PBS, pH 7.4 with 0.05% sodium azide.
Form/Appearance	Liquid
Concentration	1ug/ul
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Isotype	IgG for MAbs HP-317 and HP-387
Clonality	Monoclonal
Purity	Chromatography on protein G Sepharose
Expiration Date	12 months from date of receipt.